

doric

**Single-color Miniature Fluorescence Microscopy  
System**

User Manual

Version 3.3.2

---

# Contents

<b>1</b>	<b>Systems Overview</b>	<b>3</b>
<b>2</b>	<b>Device Overview</b>	<b>6</b>
2.1	Basic Snap-in Fluorescence Microscope Bodies . . . . .	6
2.2	Optogenetically Synchronized Fluorescence Microscope Bodies . . . . .	7
2.3	Snap-in Imaging Cannulas . . . . .	7
2.4	Dummy Microscope . . . . .	8
2.5	Pigtailed Assisted Fiber-optic & Electric Rotary Joint . . . . .	8
2.6	Fluorescence Microscope Driver . . . . .	9
2.7	Fluorescence Microscope Snapping Tool . . . . .	10
2.8	Fluorescence Microscope Holder . . . . .	10
<b>3</b>	<b>Operations Guide</b>	<b>11</b>
3.1	Connecting the Microscope Driver . . . . .	11
3.2	Installing the software . . . . .	11
3.3	Setting Up The Communication . . . . .	13
3.4	Updating The Driver Firmware . . . . .	16
3.5	Updating Doric Neuroscience Studio . . . . .	16
3.6	General Setup Guidelines . . . . .	17
3.7	Connecting the Fluorescence Microscope System . . . . .	17
<b>4</b>	<b>Using Doric Neuroscience Studio</b>	<b>24</b>
4.1	Microscope . . . . .	24
4.2	Image Analyser . . . . .	30
<b>5</b>	<b>Specifications</b>	<b>39</b>
<b>6</b>	<b>Annex 1: Cleaning and Handling</b>	<b>42</b>
6.1	Important Handling Information . . . . .	42
6.2	Cleaning Optics . . . . .	42
<b>7</b>	<b>Annex 2: Troubleshooting</b>	<b>43</b>
<b>8</b>	<b>Support</b>	<b>46</b>
8.1	Warranty . . . . .	46
8.2	Contact us . . . . .	46

## Systems Overview




Figure 1.1: Microscope System Configurations

The *Single-color miniature fluorescence microscopy systems* are configured to image a single fluorophore within small animal subject brains. An *Imaging cannula* is used to send excitation light to fluorophores and receive the resulting fluorescence. A *Miniature fluorescence microscope body* (Section 2.1, 2.2) uses a CMOS camera to take fluorescence images transmitted by the cannula and sends them to a computer. Multiple configurations and modalities are available for a multitude of specialized applications.

### 1.0.1 Basic Fluorescence Microscopy Systems




Figure 1.2: GCaMP6 Excitation Filter Set Spectra

The **Basic Fluorescence Microscopy System (BFMS)** (Fig. 1.1a) allows the imaging of a single fluorophore using a single wavelength light source. These systems are available in freely-moving and head-fixed configurations. The standard **GCaMP6** system (Fig. 1.2) uses a *Connectorized LED* as a light source, while the standard **RCaMP2** system (Fig. 1.3) uses a *Ce:YAG optical head*. Each system uses a different *Microscope Body* optimized for each fluorophore excitation/emission wavelengths.




Figure 1.3: RCaMP6 Excitation Filter Set Spectra



## 1.0.2 Optogenetically Synchronized Fluorescence Microscopy Systems




Figure 1.4: GCaMP6+NpHR3.0 Excitation Filter Set Spectra

The **Optogenetically Synchronized Fluorescence Microscopy (OSFM) system** (Fig. 1.1b) allows the imaging of a single fluorophore while performing optogenetic activation/deactivation. These systems are available in freely-moving and head-fixed configurations. When using a **GCaMP6 + NpHR3.0** microscope body a **Ce:YAG + LED optical head** is needed for illumination (Fig. 1.4). When using a **RCaMP2 + Chr2** microscope body a **Ce:YAG + Laser Diode optical head** is needed for illumination (Fig. 1.5).




Figure 1.5: RCaMP6+Chr2 Excitation Filter Set Spectra

## 1.0.3 Freely-moving and Head-fixed Configurations

The **freely-moving** (Fig. 1.1a) configuration allows fluorescence microscopy to be performed while the subject is active thanks to a *Pigtailed assisted fiber-optic & electric rotary joint*. The **head-fixed** configuration is used in combination with stereotaxic instruments. It is also used for high-precision cannula implantation.

## Device Overview

### 2.1 Basic Snap-in Fluorescence Microscope Bodies

The *SFMB Bodies* allow in vivo imaging of fluorescence within small animal subject brains. The microscopes are optimized to image calcium indicators like **GCaMP6** and **RCaMP2**. The Basic Snap-in Fluorescence Microscope Body is offered in two models: Model L for deep brain imaging and Model S for brain surface imaging. The Model L has a 0.5 NA objective lens within its body and is used with an implantable imaging cannula that transfers the image to the microscope CMOS camera. The Model S has a plan-parallel plate instead and relies on the objective lens within the model S imaging cannula to create an image on the CMOS. The *SFMB* are composed of the following elements.




Figure 2.1: Basic Snap-in Fluorescence Microscope Body

- The **Fiber-optic patch cord** connects to the **M3 optical connector**, bringing light to the microscope. The patch cord length is chosen to correspond to the electrical cable length.
- The **Electrical cable** sends and receives electrical signal for the microscope using an HDMI connector, with a standard length of 1 m. The cable is jacketed in either the light and flexible **Ultralight fiberglass jacket** or the robust but heavier **Lightweight metal jacket**.
- The **Microscope clamps** secure the microscope on an imaging cannula. This allows the removal of the microscope body between each imaging session.

- The **CMOS sensor** records fluorescence images. Each **CMOS sensor** has a serial number stored in memory that links each microscope to a specific set of mask correction filters recognizable to our software package.
- The **Imaging cannula** and **Protrusion adjustment ring** are described in Section 2.3.

## 2.2 Optogenetically Synchronized Fluorescence Microscope Bodies

The *OSFM bodies* combine fluorescence imaging and optogenetic stimulation/inhibition capabilities within the *Single-color miniature fluorescence microscope*. To avoid cross talk between optogenetic stimulation and fluorescence imaging, the *OSFM* hardware provides for at least two distinct spectral bands for optogenetic activation or fluorophore excitation (like blue and yellow). As the field of opsins and calcium indicators is very dynamic, spectral bands can be tailored to specifications. **GCaMP6 + NpHR3.0** and **RCaMP2 + ChR2** microscope versions are currently available. The appearance (Fig. 2.1) and general specifications (Table 5.1) of the *OSFM Body* are identical to the *BSFM Body*.

## 2.3 Snap-in Imaging Cannulas

Snap-in imaging cannulas transmit images of structures located inside of the brain to the surface of the skull. Cannulas have the following elements.




Figure 2.2: Snap-in imaging cannula type S (left) and type L (right)

- The **Cannula clamp groove** is used to secure the microscope body in place using **Cannula clamps**.
- The **Implantation relay lens** is found on the *Model L imaging cannula* (Fig. 2.2 right). It is an image guiding gradient-index rod lens that brings the image from inside the brain to the skull surface.
- The **Objective lens** is found on the *Model S imaging cannula* (Fig. 2.2 left). For areas near the brain surface (less than 150  $\mu\text{m}$  below the skull) it provides higher image quality and a larger FOV than the **Model L imaging cannula**.

A set of rod lenses with different lengths is available to reach different depth ranges of the brain tissue. Fine focusing of the **Model L imaging cannula** is done with the protrusion adjustment ring that comes with each cannula (Fig. 2.3). The **Model S Imaging Cannula** uses a single protrusion adjustment ring (height of 4.5 mm).




Figure 2.3: Model L Imaging Cannula Protrusion Adjustment Rings. Height from left to right: 2.05 mm, 2.77 mm, 3.48 mm, 4.20 mm and 4.92 mm

## 2.4 Dummy Microscope




Figure 2.4: *Dummy Microscope*

The **Dummy Microscope** is of similar weight and size as the microscope bodies. It is used to habituate animal subjects to the weight and feel of the microscope.

## 2.5 Pigtailed Assisted Fiber-optic & Electric Rotary Joint

To use the microscope in experiments using freely-moving animals, the *Pigtailed Assisted Fiber-optic & Electric Rotary Joint* (Fig. 2.5) is provided. This rotary joint allows effectively frictionless rotation of optical fibers and electrical cables connected to the microscope.

- The **Pigtailed Patch Cords** transmit light from the light source to the microscope. They are integrated into the rotary joint for maximal optical signal stability. The fixed and rotary patch cords use an FC connector.
- The **HDMI connectors** are used to provide electrical communication between the driver and the microscope.
- The **5 V mini-USB-B connector** connects to the power supply to allow assisted rotation.



Figure 2.5: *Pigtailed Assisted Fiber-optic & Electric Rotary Joint*

## 2.6 Fluorescence Microscope Driver

This driver allows for computer control over the excitation LED light source as well as microscope image capturing and broadcast at video rate to a single computer via high speed Ethernet communication. It can be triggered by or synchronized with external recording devices and it can also trigger other devices using 5V TTL signals (Fig. 2.6). This Fluorescence Microscope Driver is used with the *SFMB*, the *OSFM*, the *eSFMB* and the *eTFMB*.




Figure 2.6: Fluorescence Microscope Driver

- The *Fluorescence Microscope Driver* has a **12 VDC** power input (Fig. 2.6) to which the user connects the power supply provided with the system.
- The **Ethernet** port (Fig. 2.6) connects the driver to a computer using an Ethernet cable.
- The **USB-B Maintenance Port** (Fig. 2.6) is used to update driver firmware.
- The **Input BNC** (Fig. 2.6) allows the user to control the microscope using an outside, analog trigger source.
- The **Output BNC** (Fig. 2.6) provides a trigger signal corresponding to the capture of microscope images. This signal is used to trigger other devices.
- The **Microscope HDMI Connector** (Fig. 2.6) is used to connect the microscope to the driver.
- The **M8 connector** (Fig. 2.6) is used to connect the driver to a *CLED* light source.
- The **On/Off Switch** is used to turn the driver On/Off. When linking to the software, a blue circle will flash around the button. Once the driver is properly linked, the blue circle will stop flashing.

## 2.7 Fluorescence Microscope Snapping Tool




Figure 2.7: Fluorescence Microscope Snapping Tool

The Fluorescence Microscope Snapping Tool is used to remove the **Microscope clamps** from the **Cannula clamp groove**.

## 2.8 Fluorescence Microscope Holder

The Fluorescence Microscope Holder (FMH) and the Clamp for Fluorescence Microscope Holder are used to secure the microscope in a stereotaxic system.



(a) Fluorescence Microscope Holder Elements

(b) Clamp for Fluorescence Microscope Holder

Figure 2.8: Fluorescence Microscope Holder

- The *Clamp for Fluorescence Microscope Holder* (Fig. 2.8b) can be secured to a rod in a stereotaxic system using the **Slot**. The **Clamp** secures the *Fluorescence Microscope Holder* while still allowing easy removal of the holder.
- The **FMH FC Connector** allows the holder to be connected to an FC-connectorized patch cord.
- The **FMH Ferrule** is inserted into the microscope **M3 Optical Connector**. The ferrule is secured by screwing the **FMH Extremity** onto the **M3 Optical Connector**.

---

## Operations Guide

### 3.1 Connecting the Microscope Driver

- Connect the microscope driver to the power outlet using the 12 V power supply.
- Connect the microscope driver to the computer or to the router using the Ethernet cable.
- Connect the *efocus Electrical Cable* (or the HDMI cable when used with a *Pigtailed Assisted Fiber-optic & Electric Rotary Joint*) to the CMOS port on the driver.
- Push the power button. After  $\sim 5$  seconds the LED will turn ON and start to blink.

### 3.2 Installing the software

1. **Run** the Doric Neuroscience Studio Installer from the supplied USB key or download the latest version of the software from our [website](#). See Table 5.8 for computer requirements.
2. **Select** the language to use during the installation.
3. In the license agreement window (Fig. 3.1), accept the agreement and click **Next** to continue the process.




Figure 3.1: Doric Neuroscience Studio License Agreement

4. Click **Next** in the Information window.
5. **Choose** where to install the software (Fig. 3.2) and click **Next**.




Figure 3.2: Select Destination Location

6. **Choose** if desired to create a shortcut in the Start Menu folder and click **Next**.



7. **Choose** if desired to create a desktop icon and click **Next**.
8. When ready, click **Install** to begin the process. This should take a few moments. When the installation is done, the message in figure 3.3 will show up.




Figure 3.3: Successful Installation of the Doric Neuroscience Studio

9. Click **Next** and **Finish** to exit the setup.
10. Now the software is ready for use.

### 3.3 Setting Up The Communication

In order to communicate with the driver, the IP address of the computer must be static. If the driver is connected to a router, jump to section 3.3.3. If the driver is connected directly to the computer, continue to section 3.3.1.

#### 3.3.1 Configuring Static IP Address

To change the computer's IP address in Windows 7, type *network and sharing* into the Search box in the Start Menu and select Network and Sharing Center. If you are in Windows 8, it will be on the Start Screen itself (Fig. 3.4).




Figure 3.4: Open Network and Sharing Center.

On the side menu, select Change Adapter Settings (Fig. 3.5).




Figure 3.5: Click on Change Adapter Settings.

Right-click on the local adapter and select Properties (Fig. 3.6).




Figure 3.6: Right-click on local adapter Properties.

Select Internet Protocol Version 4 (TCP/IPv4) in the list, and click on Properties (Fig. 3.7).




Figure 3.7: Open IPv4 Properties.

Use the following IP address, and set the new IP address to **192.168.1.149**, and the Subnet mask to **255.255.255.0**. Leave the Default gateway and the DNS settings empty. Finally, click OK and close the Network Center (Fig. 3.8).




Figure 3.8: *Static IP settings.*

### 3.3.2 Activating The JUMBO Ethernet Frames

In order to reduce the load on the computer CPU, the Doric Neuroscience Studio is using JUMBO frames for the image transfer. In order to activate the JUMBO frames, open the Ethernet interface Properties menu as shown in Fig. 3.6 and click *Configure* (Fig. 3.9).




Figure 3.9: *Interface configuration.*

In the Configuration window, click on the *Advanced* tab, and on *jumbo frames* in the list. From the choices, select a value  $>4$  KB MTU, the bigger the better (Fig. 3.10).




Figure 3.10: Jumbo frame configuration.

### 3.3.3 Configuring The Driver IP Address

The microscope must be connected to the driver in order to proceed.

When connected for the first time, the microscope driver will broadcast its MAC address. To set the IP address, you can use the built-in ethernet setup interface in the controller software.

- Turn ON the microscope driver.
- Launch the Doric Neuroscience Studio software.
- Allow ~10 seconds for the microscope initialization.
- Choose the proper network interface from the list in the *Ethernet Settings* box, and click on **Pair** and **Remember** to configure the microscope driver.

When the system is ready to use, the power switch LED will stop blinking, and the software status will display *The microscope is ready*.

### 3.4 Updating The Driver Firmware

To update the driver firmware version, contact us at [sales@doriclenses.com](mailto:sales@doriclenses.com). The updated firmware and an installation guide will be provided.

### 3.5 Updating Doric Neuroscience Studio

To update *Doric Neuroscience Studio* see the [Neuroscience Studio Manual](#).

## 3.6 General Setup Guidelines

### 3.6.1 Optical fiber patch cord use

1. Clean the optical fiber connector before insertion. Use isopropanol and a lint-free wipe.
2. With an FC connector, the connector key must be oriented to enter within the receptacle slot to ensure proper connection (Fig. 3.11).




Figure 3.11: FC connector, Fiber Installation

**⚠** To reduce the risk of eye injury, **it is sound practice to NOT CONNECT/DISCONNECT OPTICAL FIBERS when the light source is turned on.**

### 3.6.2 Microscope Clamping Procedure

For details concerning the imaging cannula handling and implantation, as well as the microscope clamping procedure, it is important to read the [Snap-in Miniature Fluorescence Microscope Surgical Procedure](#), as well as our [Instructional Video](#).

## 3.7 Connecting the Fluorescence Microscope System

There are multiple different fluorescence microscope system configurations depending on the desired freedom of the animal and the fluorophores targeted. The animal must be held in one of the two following configuration

- The **Head-fixed** configuration (Section 3.7.1, Fig. 3.13a) allows the observation of the animal while placed in a stereotaxic apparatus.
- The **Freely-Moving** configuration (Section 3.7.2, Fig. 3.13b) allows the observation of the animal while moving freely in a cage.

Depending on the *Microscope Body* type, several different lighting configurations are used. The most common are presented in this manual.

- The **SFMB GCaMP6** configuration (Section 3.7.3, Fig. 3.14a) allows the imagery of **GCaMP6** fluorescence.
- The **SFMB RCaMP2** configuration (Section 3.7.3, Fig. 3.14b) allows the imagery of **RCaMP2** fluorescence.
- The **OSFM GCaMP6 + NpHR3.0** configuration (Section 3.7.4, Fig. 3.15a) allows the imagery of **GCaMP6** fluorescence with **NpHR3.0** activation.
- The **OSFM RCaMP2 + Chr2** configuration (Section 3.7.4, Fig. 3.15b) allows the imagery of **RCaMP2** fluorescence with **Chr2** activation.




Figure 3.12: Full RCaMP2 Freely-Moving Configuration, combination of (Fig. 3.13b) and (Fig. 3.14b)

Figure 3.13 shows the head-fixed and freely-moving system configurations. Each microscope body and lighting configuration uses the same devices in each of these configurations.

### 3.7.1 Head-fixed configuration installation

1. Install the *Fluorescence microscope holder* into the *Stereotaxic clamp* (Fig. 3.13a). Ensure the holder is tightly secured inside the clamp.
2. To secure the microscope on the *Fluorescence microscope holder*:
  - a) Remove the **Connector caps** from the microscope **M3 optical connector** and the *Fluorescence microscope holder* ferrule.
  - b) Insert the ferrule into the **M3 optical connector**. Secure them in place by screwing the *Fluorescence Microscope holder* extremity.
3. Install the *Fluorescence microscope holder* in a stereotaxic apparatus.
4. Connect the *Fluorescence microscope holder* and the light source using a fiber-optic patch cord. The patch cord characteristics (MFP\_200/230/LWMJ-0.48\_1m\_FC-FC) are set to achieve optimal coupling efficiency with the *Fluorescence microscope holder*.
5. Connect the **Pigtailed microscope cable** to the driver **HDMI connector**.
6. When ready for use, remove the cannula protective cap by unscrewing it. If using a type-L microscope, take great care to remove the cap in a slow, straight motion so as not to break the rod lens inside.




Figure 3.13: BSFM Body System Connections

### 3.7.2 Freely-moving configuration installation

1. Link the driver and the *Pigtailed assisted fiber-optic and electric rotary joint* using the HDMI cable (Fig. 3.13b).
2. Connect the top pigtailed patch cord to the light source **Beam aperture**.
3. Connect the microscope **Pigtailed HDMI Cable** to the rotary joint.
4. Connect the lower pigtailed patch cord of the rotary joint to the *FC mating adapter*, then connect the mono fiber-optic patch cord (FC-CM3) to the microscope. The patch cord (MFP\_200/230/LWMJ-0.48\_0.8m\_FC-CM3) linking the Pigtailed Assisted Fiber-optic & Electric Rotary Joint (AHRJ) and the microscope is chosen to ensure proper illumination characteristics.
5. Remove the **Input Protective Cap** from the cannula and insert the microscope body.

### 3.7.3 Connecting the BFMS

The *BSFM body* comes with two standard models; the first allows **GCaMP6** fluorescence imagery while the second allows **RCaMP2** fluorescence imagery.




Figure 3.14: BFMS Body Illumination System Connections



## Microscope Driver and Light Source Installation


**⚠** If using the *Ce:YAG Optical head*, connect the **Interlock connector plug** to a **Laser safety interlock circuit** for safe use. See the application note [Important Laser Safety Information](#) for more information, or contact your institution's laser safety officer.

1. Connect the *12 VDC power supply* to the microscope driver **12 VDC power input**.
2. Connect the driver and the light source.
  - If using an LED light source (such as for the **GCaMP6** system (Fig. 3.14a)), connect the driver to the LED using the **M8 connector**.
  - If using a *Ce:YAG Optical Head* (such as for the **RCaMP2** system (Fig. 3.14b)), see the **Ce:YAG manual** for details on driver and light source installation. Once complete, connect the *Fluorescence microscope driver BNC Output* to the *Ce:YAG Driver BNC Input*. The *Ce:YAG Driver* must be put in **TTL** mode to receive proper signals from the *Microscope driver*.
3. Connect the driver to the computer using the *CAT5E Ethernet Cable*.
4. Turn on the *Fluorescence Microscope Driver*.


### 3.7.4 Connecting the OSFM system

The indications in section 3.7, notably Fig. 3.13b and Fig. 3.13a are to be followed for the installation of the **OSFM**. As the **OSFM** uses two light sources, driver and light source installation is more elaborate. Both standard configurations use the *Ce:YAG + LED/LD fiber Light Source* to illuminate the sample<sup>1</sup>.

#### Driver and light source installation



(a) GCaMP6 + NpHR3.0 configuration



(b) RCaMP2 + ChR2 configuration

Figure 3.15: OSFM microscope system connections

<sup>1</sup>For details and specifications concerning the Ce:YAG + LED/LD Fiber Light Source, see their manual.

### **GCaMP6 + NpHR3.0 configuration**

For the **GCaMP6 + NpHR3.0 configuration** (Fig. 3.15a), the *Ce:YAG + LED optical head* is used as the light source. The LED light source is used for excitation, while the *Ce:YAG* light source is used for opsin activation.

1. Connect the *OTPG* to the computer using a **USB-A/USB-B cable**.
2. In this configuration, the LED light source needs to be triggered by the *Fluorescence microscope driver*. Connect the driver **BNC output** to the *Ce:YAG + LED driver* **LED BNC input**.
3. Connect the *Ce:YAG + LED driver* **Ce:YAG BNC input** to an *OTPG* channel.
4. Connect the *Fluorescence microscope driver* **Input BNC** to an *OTPG* channel.
5. When installing in a head-fixed or freely-moving configuration, the rest of the system is connected as described in section 3.7.1 and 3.7.2.

### **RCaMP2 + ChR2 configuration**

For the **RCaMP2 + ChR2 configuration** (Fig. 3.15b), the *Ce:YAG + LD optical head* is used as the light source. The LD light source is used for opsin activation, while the *Ce:YAG* light source is used for excitation.

1. Connect the *OTPG* to the computer using a **USB-A/USB-B cable**.
2. In this configuration, the *Ce:YAG* light source needs to be triggered by the *Fluorescence microscope driver*. Connect the driver **BNC output** to the *Ce:YAG + LD driver* **Ce:YAG BNC input**.
3. Connect the *Ce:YAG + LD driver* **LD BNC input** to an *OTPG* channel.
4. Connect the *Fluorescence microscope driver* **Input BNC** to an *OTPG* channel.
5. When installing in a head-fixed or freely-moving configuration, the rest of the system is connected as described in section 3.7.1 and 3.7.2.

## Using Doric Neuroscience Studio

### 4.1 Microscope

The Microscope module of the Doric Neuroscience Studio provides an interface to control our Fluorescence Microscope Driver. The module enables image acquisition and its export in 16 bit .tif or in .doric (hdf5-based) files. The TIF format can easily be read with any standard imaging software. Doric files can be read by the Doric Neuroscience Studio **Image Analysis Module** or using an HDF5 library. Despite the fact that the images are saved with a 16 bit pixel depth, the true image pixel depth is 10 bit, so pixel gray values are contained between 0 and 1020 counts.

Below is the user interface (Fig. 4.1) and a complete description of all the functions.




Figure 4.1: Microscope Module Interface

1. The **Image Box** (Fig. 4.1) displays images from the microscope and allows region of interest (ROI) drawing by clicking and dragging the mouse over the image.

- The **Sensor Tabs** (Fig. 4.1) display the sensors available to view. For multi-sensor microscopes, changing tabs allows you to see the image available to each.
- The microscope **Status** (Fig. 4.1) will indicate the current microscope state (Live/Stopped).
- The **Exposure** (in ms) (Fig. 4.1) indicates the exposure time of the microscope sensor.
- The **FPS** (Frames Per Second) (Fig. 4.1) indicates the number of frames per second taken by the sensor.
- The **Gain** (Fig. 4.1) indicates the electrical gain of the sensor.
- The **Size** (Fig. 4.1) indicates the resolution of the sensor images (in Pixels x Pixels).
- The **Bin** (Fig. 4.1) status indicates whether or not the sensor image is being binned (yes/no).

2. The **Capture** tab (Fig. 4.1) contains different image-capturing functions of the microscope.



Figure 4.2: Capture Tab

- The **Live** button (Fig. 4.2), when pressed, displays images from the microscope. These images are not saved.
- The **Snap** button (Fig. 4.2), when pressed, takes a snapshot of the current image and saves it in the requested directory with the desired name (**Saving Options**) as a single image.
- The **Album** button (Fig. 4.2), when pressed, acquires a snapshot and adds it to an album stack. The whole stack can be saved as one image stack.
- The **Record** button (Fig. 4.2), when pressed, acquires a continuous image stream, until **Stop** is pressed, and saves it in the requested directory with the desired name (**Saving Options**) as one image stack.




Figure 4.3: Time Series Window

- The **Time Series** button (Fig. 4.2), when pressed, opens the time series interface (Fig. 4.3).
  - The **Number of time points** (Fig. 4.3) defines the number of moments when a set of images will be recorded.
  - The **Images per time point** (Fig. 4.3) defines the number of images taken in each set.
  - The **Time interval between points** (Fig. 4.3), defined in ms, s and min, defines the duration between each image set. This duration always has a minimum value of **Exposure time x Images per time point**.
  - The **Summary** box (Fig. 4.3) shows many values related to the time series, including the **Total images** recorded, the **Total memory** occupied by the full series, the **Time point duration** and the **Total duration** of the full series.
  - The **Progression bar** (Fig. 4.3) displays the progress (in %) of the time series.





Figure 4.4: Saving Options Window

- f) The **Saving options** button opens the **Saving options window**.
  - i. The **Filename** box (Fig. 4.4) is used to define the recorded file name.
  - ii. The **...** button (Fig. 4.4) opens a window used to choose the save file location.
  - iii. The **File type** drop-down menu (Fig. 4.4) is used to decide which file type is used to save images. For files larger than 4 GB, the .doric extension is recommended.
  - iv. The **Index** box (Fig. 4.4) displays the current index that will be added to the filename.
  - v. The **Target File** box (Fig. 4.4) shows the full location and name of the file being saved when an image sequence is recorded.

3. **Microscope settings** tab (Fig. 4.1) is used to set parameters related to the microscope recording images.

[Microscope Settings Tab (2-Color Fluorescence Microscope)]



[Microscope Settings Tab (eFocus Miniature Fluorescence Microscope)]






Figure 4.5: Microscope Settings Tab


- a) The **Exposure** box (Fig. 4.5) sets the exposure time of the sensor. The time can be set between 22 and 1000 ms.
- b) The **SENSOR** section (Fig. 4.5) defines characteristics for a single sensor and the associated excitation source. When a microscope used has multiple sensors, multiple **SENSOR** sections will be displayed, one for each sensor.
  - i. The **Gain** box defines the sensor gain.
  - ii. The **Illuminator power (%)** box defines the power emitted by the excitation light source. The light sources will be activated when the image acquisition is started. The maximum optical power (in mW) depends on the light source model.
  - iii. The **Illumination source** box defines if a light source is linked to a given sensor. If No source is selected, this light source can be controlled independently using the light source Tab. For more information on how to control the light sources using the Light Source Tab, please refer to the [Doric Neuroscience Studio](#) User Manual.
- c) The **Working Distance Adjustment Slider** appears when an eFocus Miniature Fluorescence Microscope is connected to the driver. This slider will adjust the working distance from -45 to 45  $\mu\text{m}$  for snap-in fluorescence microscope bodies and from 0 to 350  $\mu\text{m}$  for twist-on fluorescence microscope bodies.



(a) External Trigger Settings Window



(b) Gated Mode



(c) Trigger Out Mode

Figure 4.6: External Trigger Settings Windows

- d) The **External Trigger** button (Fig. 4.5) opens the **external trigger** window.
- i. The **Number of images per trig** box (Fig. 4.6a) defines the number of images acquired at each trigger pulse.
  - ii. The **File name/location** (Fig. 4.6a) box displays the location where the images are saved as well as their file name.
  - iii. The **Select...** (Fig. 4.6a) button allows the selection of the **File name/location**.
  - iv. The **Progression bar** (Fig. 4.6a) displays the advancement of the triggered sequence (in %).
  - v. The **Gated mode** checkbox (Fig. 4.6a) will change the external trigger to gated mode (Fig. 4.6b). In this mode, the microscope will only acquire images when a high TTL signal is received on the TRIG IN input.
  - vi. Selecting the **Trigger Out** Tab will change the external trigger to *Trigger Out* mode. This mode is used to select the type of TTL signal generated at the output of the Microscope Trigger Out BNC. In Normal Mode, the TTL signal is high whenever the microscope is live and in Triggered with each frame Mode, a TTL pulse is generated at each frame.
- e) The **Save configuration** button (Fig. 4.5) will save all **Microscope settings** and **Image settings** in a **.doric** format file.
- f) The **Load configuration** button (Fig. 4.5) will load a selected configuration file.
- g) The **Select mask file** button opens a window to select a mask file for the microscope used. This section only appears when a *2-color Fluorescence Microscope* or an *efocus Microscope* is connected. The mask file currently loaded will be shown just above it. For more information on masks, see section ??.




Figure 4.7: Image Settings Tab

4. The **Image settings** tab (Fig. 4.7) is used to define certain settings related to the displayed and recorded images.
  - a) The **Crop Image** button (Fig. 4.7) allows a square to be drawn onto the image. When a new **Capture** sequence is activated, only the cropped region will be captured.
  - b) The **Reset crop** button (Fig. 4.7) resets the cropped image to its original state. The change will only appear when a new **Capture** sequence is activated.
  - c) The **Binning** drop-down list (Fig. 4.7) allows the binning of pixels. This reduces the number of pixels for smaller save file sizes.
5. The **View** tab (Fig. 4.1) is used to change viewing parameters of the sensor image. These changes will only appear on the sensor image when a new **Capture** sequence is started. Any adjustments made affect only the displayed image and not the recorded images.




Figure 4.8: View Tab

- a) The **Zoom In/Zoom Out** buttons (Fig. 4.8) will increase/decrease the zoom of the sensor image.
- b) The **Reset Zoom** button (Fig. 4.8) will reset the **Zoom factor** to 100%.
- c) The **Zoom Factor** drop-down list (Fig. 4.8) allows the selection of a zoom factor from a pre-set list. The box will also display the current zoom if it was changed using different buttons.
- d) The **Fit In View** button (Fig. 4.8) resize the image to fit the size of the Image Viewer box.
- e) The **Roi shape** drop-down list (Fig. 4.8) allows the selection of the shape used when drawing a **Region Of Interest** onto a sensor image. These shapes include **Freehand, Circle, Rectangle** and **Square**.
- f) The **SENSOR** section (Fig. 4.8) is used to adjust contrast on a given sensor. When a microscope used has multiple sensors, multiple **SENSOR** sections will be displayed, one for each sensor.
  - i. The **Contrast** slider (Fig. 4.8) allows the adjustment of contrast from 0.1 to 5.
  - ii. The **Min/Max** sliders (Fig. 4.8) indicate the minimum/maximum number of counts displayed. Should the **Min** be above 0, all pixels with lower count will display a minimal value. Should the **Max** be below 1020, all pixels with a higher count will appear saturated.
  - iii. The **Auto contrast** slider button (Fig. 4.8) will active an automatic contrast adjustment algorithm.
  - iv. The **Reset** button resets contrast functions to their default settings.
- g) The **Pseudocolor** drop-down list (Fig. 4.8) allow the sensor image color palette to be changed.
- h) The **Show saturation** checkbox (Fig. 4.8) allows all saturation on the sensor image to be displayed in red. This function is only available if no pseudocolor is selected.






Figure 4.9: Microscope Ethernet Tab

6. The **Ethernet** tab (Fig. 4.1) is used to define the ethernet connection used to connect the computer to the microscope driver.
  - a) The **Refresh** button (Fig. 4.9) will identify any accessible IP addresses and add them to the drop-down list.
  - b) The **Ethernet** drop-down list (Fig. 4.9) includes all IP addresses connected to an ethernet adapter. The proper one must be selected to properly connect the microscope.
  - c) The **Pair** button (Fig. 4.9) connects the software to the driver.
  - d) The **Remember** checkbox (Fig. 4.9) will keep the chosen IP address so that the chosen microscope driver will be connected automatically next time the software is opened.



Figure 4.10: ROI Manager

7. The **ROI Manager** (Fig. 4.1) displays the live mean pixel intensity from a drawn ROI.
  - a) The **Average Intensity in ROI** plot (Fig. 4.10) displays the average intensity over time inside a drawn ROI. *CTRL + mouse wheel* will adjust the x-axis zoom, while *SHIFT + mouse wheel* will adjust the y-axis zoom.
  - b) The **Clear** button (Fig. 4.10) will clear any data displayed in the ROI manager and the ROI on the **Image Viewer**.
  - c) The **Optimal Zoom** button (Fig. 4.10) sets the zoom factor on the plot to best display all data.
  - d) The **Reset Zoom** button (Fig. 4.10) resets the zoom to its default setting.

To update the driver firmware version, connect the driver to the computer using a USB-A/USB-B male cable **only**. In the **Hardware** drop-down menu, select the **Update Microscope** to start the process. Turn the driver power ON.

## 4.2 Image Analyser

This module provides an easy way to extract relevant data from the images acquired by the Doric miniature fluorescence microscopes. The software loads images in .TIF and .doric formats, implements image processing functions and an export tool to save the fluorescence data in .CSV or .doric format. This software does not replace standard analysis tools such as Matlab, ImageJ or Excel, but aims to offer useful processing algorithms developed for the microscope images. All the underlying algorithms are implemented from the [OpenCV](#) library. In this section, we will describe the different functions available, and how to use them.




Figure 4.11: Image Analysis Module Interface

1. The **Image Viewer** displays the loaded images, allows navigation through the image stack and the drawing of regions of interest (ROIs) by clicking and dragging the mouse over the image. Multiple image sets can be opened, appearing as tabs in the upper left of the image box.
2. The **ROI Manager** displays the different ROI parameters and traces the mean signal intensity over time for each ROI.
3. The **Function Toolbar** contains all the buttons and functions accessible.

### 4.2.1 Function Toolbar

1. The **File tab** (Fig. 4.12) is used to save/load data.




Figure 4.12: File Tab

- The **Load Images** function loads a square, 16 bit .tif file or .doric file.
- The **Save Images** function saves the current image tab to a 16 bit TIF multipage file or .doric file.
- The **Export Traces** function saves the average fluorescence intensity values for each ROI of the current tab to a .CSV or .doric file.
- The **Export Raster Image** function saves the raster plot of the ROIs of the current image tab to a .jpg file.

2. The **ROI tab** (Fig. 4.13) is used to save/load data relating to regions of interest drawn on an image.




Figure 4.13: ROI Tab

- The **Load ROIs** function loads .CSV file containing informations about the saved ROIs.
- The **Save ROIs** function saves the current ROIs information to a .CSV file.
- The **Clear All** button clears all ROIs.
- The **ROI shape** function is a drop-down list that allows the selection of the **ROI** shape. These include **Freehand, Circle, Rectangle** and **Square**.

3. The **Spikes tab** (Fig. 4.14) is used to save/load data relating to regions of interest drawn on an image.




Figure 4.14: Spikes Tab

- The **Load Spikes** function loads .CSV file containing informations about the saved Spikes.
- The **Save Spikes** function saves the current Spikes information to a .CSV file.
- The **Clear All** button clears all Spikes.

4. The **View tab** (Fig. 4.15) is used to manipulate the appearance of an image without changing base data.




Figure 4.15: View Tab

- The **Fit In View** check box adjusts automatically the size of the current image to the **Image Viewer** window.
- The **Reset Zoom** and **Zoom factor** functions adjusts the display size of the current image.
- The **Contrast** function applies a different luminance response curve (gamma). See section 4.2.3 for details.
- The **Min** function applies a lower threshold with the cut-off value defined by the slider. See section 4.2.3 for details.
- The **Max** function applies an upper threshold with the cut-off value defined by the slider. See section 4.2.3 for details.

- The **Autocontrast** function directly applies the `equalizeHist` function of the OpenCV library.
- The **Reset** function returns the contrast and range values to their default.
- The **Image Info** button displays the image information window.
- The **Pseudocolor** function is a drop-down list for selecting alternate coloring schemes for the images presented.
- The **Frame Display Time** function adjust the frame rate in **Play** mode.

5. The **Processing tab** (Fig. 4.16) is used to process the image data.




Figure 4.16: *Processing Tab*

- The **Crop** function allows to crop the current image in smaller dimensions to reduce the amount of data and facilitate the processing.
- The **Bin** function combines a cluster of pixels in a single pixel to reduce the amount of data and facilitate the processing. Note: in 2x2 binning, an array of 4 pixels becomes a single larger pixel.
- The **Discard Frame** function allows to remove user-defined frames in a data set. Note: The timestamps of the remaining frames stay the same when discarding frames.
- The **Align** function aligns the image stack to the user-defined key frame. See section 4.2.3 for computational details. Selecting this button will open the **Align Images** window (Fig. 4.17). By selecting the **Save Alignment Values** checkbox, the image alignment values will be preserved when saving the processed images. There are 4 different methods available.





Figure 4.17: *Align Images Window*

- The **First Frame Of Current Image Set** method uses the first image in the set to align the rest.
- The **Select Frame From Current Image Set** method allows the selection of a single image in the set to use for alignment of all other frames.
- The **Select Other Image Set And Frame** method aligns the current set using data from a different image set.
- The **Select From Alignment Value File** method uses a previously-defined alignment for another image set. This method is most valuable when trying to align images from the *2-color fluorescence microscope*, to align one color channel using the data from the other.
- The **Remove Background** function removes the average value of a selected ROI from all images in the stack. Note: it is not recommended to use the **Remove Background** function before with the  $\Delta F/F_0$  function.
- The  $\Delta F/F_0$  function calculates the normalized fluorescence variation of the images and displays the results in a new tab. When selected, the See section 4.2.3 for details.
- The **Find Cells** function detects the cells and creates the ROI automatically. See section 4.2.3 for details.
- The **Stack Projection** function projects all movie frames to a single frame using the method selected in the Settings dialog. See section 4.2.3 for details.

- The **Find Spikes** function detects the spikes of the traces calculated from the ROI. The positions of the spikes are indicated by a red dot in the **ROI Manager**. See section 4.2.3 for details.
- The **Batch Processing** function opens the **Batch Processing Window** (Fig. 4.18). This allows the processing of large datasets in sequential order, without needing to activate each individual function. The processing defined in the batch processing window is applied to all the data saved in the destination file.



(a) Batch processing window

(b) Typical batch processing sequence

Figure 4.18: Batch Processing Window

- a) The **Available processes** box lists all processes available. Processes on the list will be greyed out if the work-flow order prevents them from being used. Each process has a number of parameters that are identical to those used outside of batch processing.
- The **Crop** function allows to crop the current image in smaller dimensions to reduce the amount of data and facilitate the processing.
  - The **Bin** function combines a cluster of pixels in a single pixel to reduce the amount of data and facilitate the processing.
  - The **Convert** process is used to convert an image stack to **.doric** or **.tif** format.
  - The **Align** process aligns the image stack to the user-defined key frame. See section 4.2.3 for computational details.
  - The  **$\Delta F/F_0$**  process calculates the normalized fluorescence variation of the images and displays the results in a new tab. See section 4.2.3 for details.

- The **Find Cells** process detects the cells and creates the ROI automatically. See section 4.2.3 for details.
  - The **Stack Projection** process projects all image frames to a single frame using the method selected in the Settings dialog. See section 4.2.3 for details.
- b) The **Workflow** box displays the order in which image processing actions will be taken. The parameters of the selected functions are adjusted in the **Options** box.
  - c) The **Choose a Folder** button allows the selection of a folder to save batch processing results.
  - d) The **File Type** list is used to defined the file extension used when the images are saved.
  - e) The **Save intermediate files** option will save intermediary files in the image processing process alongside the completed files.

## 4.2.2 ROI Manager





Figure 4.19: ROI Manager


The ROI manager extracts average intensity of a defined section of the image over an entire image stack. There is no limit to the number of ROI allowed per image stack.

1. The **Image Viewer** contains the image stack and the ROI, numbered according to the order they were set. The ROI can be saved independently from the image stack on the ROI toolbar. The ROI are drawn directly on the *Image Viewer* in a *freehand* manner. All selected ROI can be moved together directly in the *Image Viewer*.
2. The **Intensity Plot** panel shows the plot of average intensity as a function of the frame index. The Y-axis represents the average count of all the pixels of the ROI. It is separated into *Overview* and *Detailed* tabs.
  - The **Overview** tab displays all the traces on the same graph, on the same scale (see Fig. 4.20a).
  - The **Detailed** tab displays each trace on a separate graph, allowing for precise intensity measurements (see Fig. 4.20b).
  - The **Raster** tab displays all the traces on the same color-coded graph. In a raster plot each row (y-axis) corresponds to an ROI. The columns (x-axis) correspond to the current time (see Fig. 4.20c).
3. The **ROI Data** list shows the parameters defining each ROI. Selected items will be displayed in orange on the *Image Viewer* and in the *Overview* graph.
  - The **ID** shows the order of the ROI (starting at 0).
  - The **Area** shows the area (in pixels) contained in the ROI.






(a) Overview Graph



(b) Detailed Graph



(c) Raster plot

Figure 4.20: ROI Manager Graph Tabs

### 4.2.3 Algorithms

#### Contrast

The contrast adjustment applies the following operation to each pixel of the image:  $V_{out} = AV_{in}^\gamma$ , where  $V_{out}$  is the corrected pixel value,  $A = 1$ ,  $V_{in}$  is the initial pixel value, and  $\gamma$  is the value as selected by the contrast slider.

#### Min and Max ranges

When the values of the display range are other than the default  $min = 0$  and  $max = 1020$ , the following operation is applied to each pixel:  $V_{out} = 1020 * (V_{in} - min) / (max - min)$ , where  $V_{out}$  is the corrected pixel value,  $V_{in}$  is the initial pixel value,  $min$  and  $max$  are respectively the minimum and maximum slider values.

#### Image Alignment

The algorithm is inspired from Manuel Guizar-Sicairos, Samuel T. Thurman, and James R. Fienup, *Efficient subpixel image registration algorithms*, Opt. Lett. 33, 156-158 (2008). The basic idea is to obtain an initial estimate of the crosscorrelation peak by a Fourier transform and then refine the shift estimation by upsampling the Fourier transform only in a small neighborhood of that estimate by means of a matrix-multiply Fourier transform. With this procedure, all the image points are used to compute the upsampled crosscorrelation. In order to increase the precision of the algorithm, we use the laplacian of the images as inputs, instead of using the raw images. Briefly, the algorithm applies the following steps:

1. Calculate gaussian blur of the reference image with window of size 39 to smooth high frequency noise.
2. Calculate the laplacian of the blurred reference image.
3. Use the absolute values as the final reference image.

4. Reproduce steps 1 to 4 for the following image.
5. Calculate the 2D Fourier transform of the reference and the target image.
6. Multiply both images.
7. Calculate the inverse Fourier transform of the product image.
8. Get the position of the maximum correlation peak.
9. Create an upsample array around the maximum correlation peak to refine the shift calculations.
10. Calculate the Fourier transform of the larger array.
11. Do the matrix multiplication.
12. Locate the maximum correlation and map it back to the original space.

### $\Delta F/F_0$

The algorithm calculates a standard  $\Delta F/F_0$  with  $F_0$  corresponding to the temporal average intensity, with an optional preprocessing step to remove the illumination variation artefacts. In order to properly calculate the  $\Delta F/F_0$ , the algorithm uses a dark frame to account for the sensor electronic offset. Calculating the  $\Delta F/F_0$  without subtracting the offset will lead to artificially lower values. To record a dark frame, set the microscope driver to the desired exposure and gain, the LED power to zero and take a snapshot. Before calculating the  $F_0$ , the average temporal variations can be compensated to get a flat temporal average profile (Fig. 4.21). Keep in mind that removing the average temporal profile can also remove global activity patterns.




Figure 4.21:  $\Delta F/F_0$  Settings

Briefly, the algorithm applies the following steps:

1. Calculate the average image intensity as a function of time ( $C$ ).
2. If the global variation removal option is selected, apply the following correction to each image:  $I_{out} = (I_{in} - I_{dark}) * (mean(C - I_{dark}) / (C - I_{dark}))$  where  $I_{out}$  is the LED illumination corrected image,  $I_{in}$  the input image and  $C$  is the average temporal trace.
3. Calculate  $F_0$  as the average projection of the movie.
4. Calculate the relative change  $R(t)$  of fluorescence signal  $R(t) = (F(t) - F_0) / F_0$ .

### Find Cells

The algorithm is inspired by Eran A. Mukamel, Axel Nimmerjahn and Mark J. Schnitzer, *Automated analysis of cellular signals from large-scale calcium imaging data*, Neuron 63(6), 747-760 (2009). The basic idea is to use a principal component analysis (PCA) as input of an independent component analysis (ICA) to separate the different temporal signals contained in the movie. This method is used as a starting point to determine the position of the different active cells. It is coupled with a segmentation routine optimized for reducing the false positives. The *Find Cells* algorithm uses user-defined boundaries shown in Fig. 4.22. The first parameter is an estimate of the number of cell present in the movie. By design, it must be lower than the number of frames minus five. The next parameters are the smallest and biggest object diameter in microns. These values are used to filtered the object found by the PCA/ICA.






Figure 4.22: *Find Cells Settings*

Briefly, the algorithm applies the following steps:

1. Calculate and remove the spatiotemporal average from the movie, as the PCA/ICA algorithm requires zero-mean data.
2. Run OpenCV PCA algorithm on the centered data.
3. Normalize data by standard variation.
4. Calculate ICA with PCA as input data.
5. Apply segmentation to each ICA found.
6. Filter contours found at the previous step using user-defined boundaries.

### Stack Projection

This function can be used to help for ROI drawing. It calculates a temporal projection using the user-defined method (see Fig. 4.23).




Figure 4.23: *Stack Projection Settings*

**Maximum:** the output is the maximum value found in all frames for each pixel.

**Average:** the output is the mean value of all frames for each pixel.

**Sum:** the output is the sum of all frames for each pixel.

**Minimum:** the output is the minimum value found in all frames for each pixel.

## Find Spikes

This function detects the spikes of the traces calculated from the ROIs. The detection threshold is defined as a factor of the standard deviation of the trace.




Figure 4.24: *Spike finder Settings*

## Specifications

Table 5.1: *Single-color Fluorescence Microscope Body Specifications*

SPECIFICATION	Microscope Body	
	Model S	Model L
Mass without cables (g)	2.2	
Dimensions without cables in mm (W x L x H)	8.8 x 13.9 x 16.6	
Frame rate (fps)	49	
Objective lens NA	0.5	
FOV at image plane (pixel)	630 x 630	
FOV at object plane ( $\mu\text{m}$ )	700 x 700	350 x 350
Lens magnification	3.3x	6x

Table 5.2: *Single-color Fluorescence Microscope Body Excitation and Detection Spectra*

Microscope Bodies	SPECTRUM (nm) <sup>1</sup>		
	Excitation	Detection	Opsin activation
<b>Basic</b>			
GCaMP6	445/62	525/40	—
RCaMP2	549/15	609/57	—
<b>OSFM</b>			
GCaMP6 + NpHR3.0	445/62	525/40	616/76
RCaMP2 + ChR2	549/15	609/57	Compatible with 450, 473, 488

Table 5.3: *Model S Cannula general specifications*

SPECIFICATIONS	VALUE
Magnification	2 x
Working distance	2.4 mm

<sup>1</sup>Center wavelength/bandwidth

Table 5.4: Model L Cannula general specifications

SPECIFICATIONS	VALUE
Lens diameter	500 $\mu\text{m}$
Working distance	80 $\mu\text{m}$

Table 5.5: Fluorescence Microscope Driver general specifications

SPECIFICATIONS	VALUE	NOTES
Power supply	110 - 240 VAC, 50 - 60 Hz	
DC Power supply	12 VDC	
Dimensions	186 x 90 x 77 mm <sup>3</sup>	Including connectors
Data link	Gigabit ethernet	
TTL input voltage	0 to +5 V	
LED Analog input voltage	400 mA/V light source current 40 mA/V light source current	Standard 1 A LED Driver Low power mode enabled
LED BNC output voltage	2.5 V/A	
LED Maximum output current range	200, 2000 mA	
LED Maximum forward voltage	7 V	
LED Minimum output current	2.5 mA	Low power mode enabled
LED Rise/Fall time	<10 $\mu\text{s}$	Typical
LED Connector	M8	See Figure 5.1



Figure 5.1: M8 Female Pinout (Microscope Driver)

Table 5.6: Fluorescence Microscope Driver Software specifications

SPECIFICATIONS	VALUE	NOTES
Current adjustment steps	1 mA	
Modulation minimum frequency	0.01 Hz <sup>2</sup>	Internal complex mode : 0.000054 Hz
Modulation maximum frequency	50 kHz	-3 dB attenuation
Minimum ON or OFF time	0.005 ms <sup>2</sup>	Internal complex mode : 2 ms
Maximum ON or OFF time	100 s <sup>2</sup>	Internal complex mode : 5 h
Maximum number of pulses per sequence	16.68 millions <sup>2</sup>	Internal complex mode : 65 535
Maximum number of sequences	4.2 billions <sup>2</sup>	Internal complex mode : 65 535
Minimum step increments	39 $\mu\text{sec}$ <sup>2</sup>	Internal complex mode only
Number of steps per period	128 <sup>2</sup>	Internal complex mode only
Scope Acquisition speed	10 kS/s	Single channel

<sup>2</sup>For all operation modes, except the internal complex mode

Table 5.7: Connectorized LEDs general specifications

SPECIFICATION	VALUE	NOTE
Input current	Min: 0 mA - Max: 1000 mA	Recommended 700 mA
Forward voltage	3.0 - 4.0 V Typical	
Dimensions	63 x 59 x 21 mm <sup>3</sup>	Without cable, with baseplate
Mass	~100 g	
CLED cable connector	M8-4pins - Male	See Figure 5.2
Output NA	0.55	NA of up to 0.63 NA will slightly increase power;
Output optical fiber core diameter	<960 μm	Power scales up to this core diameter




Figure 5.2: M8 Male Pinout (CLED)

Table 5.8: Doric Neuroscience Studio Hardware Requirements

SPECIFICATIONS	VALUE	NOTES
Operating System	Windows 7, 8, 10	64-bit
Memory (Minimum/Recommended)	4 GB/16 GB	
Processor Speed (Minimum/Recommended)	2 Ghz Quad-Core i5/ 3.46 Ghz Eight-core i7	
Hard Drive	500 MB	

## Annex 1: Cleaning and Handling

### 6.1 Important Handling Information

*Warning:* Handle the microscope and the cannula with care.

Miniature fluorescence microscopes are composed of sensitive electronic components and should always be handled with care. When they are not used, the microscope body with its cannula should be stored in a closed, dust-free environment. Some microscope components must be handled with extra care:

- **Electrical cable: Do not twist or pull on the cable.** This cable is pigtailed to the CMOS sensor and cannot be easily replaced.
- **Relay lens:** The cannula lens is made of glass and is unprotected. **Abrasive materials can scratch the surface** and reduce the image quality.

The microscope bodies and the cannula lenses are made of glass, metal, plastic and the contact with organic tissues or liquids, like blood or salted water solution is not recommended. If the microscope comes in contact with these substances, clean the optics (section 6.2) to avoid the apparition of stains.

### 6.2 Cleaning Optics

The microscope objective lens should be cleaned before each use. The procedure explained here can also be used to clean the cannula relay lenses.

- Turn the driver OFF.
- **Wear gloves to manipulate the microscope.** Finger oil can stain the glass and is often hard to remove properly.
- Use isopropyl alcohol on a cotton swab to gently clean the lens.
- **Do not blow on the optics.** Saliva particles will often stain the surface. Larger dust particles can be removed using a dust-free blower before cleaning with a cotton swab.

## Annex 2: Troubleshooting

### 7.0.1 Software

How to make the *Doric Neuroscience Studio* detect the *Microscope Driver*?

1. Ensure that the *Microscope Driver* is plugged into the computer using an Ethernet cable.
2. Ensure that the each *Electrical Cable Connector* is plugged into the appropriate device. The *Microscope Driver* must be linked to the *Microscope*.
3. Ensure that the IP address is static (see section 3.3.1)
4. Ensure the *Jumbo Frames* are activated (see section 3.3.1)
5. *Windows Firewall* can prevent communication. To ensure the communication is not being blocked, open the *Windows Firewall* configuration window, then click on *Allow an app through the firewall*. From there, select the *Change Settings* button, find the *Doric Neuroscience Studio* and check the *Private* and *Public* checkboxes.
6. In the *Network & Sharing Center*, check the Ethernet connection; it should indicated *Unidentified Network*. If *Network Cable Unplugged* is shown despite the Ethernet cable being plugged in and the driver being turned on, disable and re-enable the Ethernet connection.
7. Ensure *Network & sharing* is properly configured at 1 Gbps by double-clicking the Ethernet connection and checking the *Speed*.
8. When the *1-color Microscope Driver* is activated, the *On/Off Switch* should blink blue while initializing. If the light is sustained without any blinking when first turned on, restart the *Microscope Driver*.
9. Certain Intel Ethernet cards must be activated in *Slave Mode* to function. This mode can be found in the same menu as the *Jumbo Frames* (see section 3.3.2).

How can I stop the software lagging and/or dropping frames<sup>1</sup>?

1. Deactivate all internet using programs that can conflict with the *Doric Neuroscience Studio* (IE Skype, Firewall, etc.)
2. Use a computer with the recommended specifications:
  - *Operating System*: Windows 10
  - *CPU*: Quad Core I7 3.46 GHz
  - *RAM*: 16 Gb

---

<sup>1</sup>Dropped frames are black frames that occur when an image is lost in communication. They can easily be spotted in the *Average Intensity In ROI* trace if the value descends to 0.

- *Dedicated Graphics Card*: with Open GL version 4.6 recommended
  - *Desktop computer recommended*
3. Windows might limit the Ethernet performances to reduce energy consumption. To ensure that the communication is not limited, open the Power option window:
- Press the Windows + R keys to open the Run dialog box.
  - Type in the following text: "powercfg.cpl", and then press Enter.
  - In the Power Options window, under Select a power plan, choose High Performance.
  - If you do not see the High Performance option, click the down arrow next to Show additional plans.
  - If available, change the System standby and System hibernates settings to Never.
  - Click Save changes or click OK.

## How can I visualize recorded frames?

1. All images will appear black in *Window Image Preview*/Traditional Image Viewers as they are special 16 bit .tif files. Use dedicated software such as the *Doric Neuroscience Studio Image Analyzer* or *ImageJ*.
2. Due to library usage conflicts, the *Dell Backup & Recovery* application interferes with the loading of images in the *Doric Neuroscience Studio*. Uninstall the *Dell Backup & Recovery* application from the computer

## Can I use a USB to Ethernet adapter to connect the driver?

1. The microscope driver must be connected to a computer ethernet port.
2. Should a USB to Ethernet adapter be used for other function, such as internet access, the adapter must be disabled during the first initialization of the microscope.

### 7.0.2 Hardware

## How do I prevent instability in the *Assisted Opto-electric Rotary Joint*?

1. Ensure the optical fiber *Patch Cord* is of equal or shorter length than the microscope *Electrical Cable* when connector to the *Assisted Opto-Electric Rotary Joint*. Even if the cable is looped, the distance from rotary joint to patch-cord connector should be shorter than the length of the electrical cable.

## How do I prevent the *Cannula* turning in the *Protrusion Adjustment Ring*?

1. These two components are meant to be glued together after installation. If they have not been glued during installation, add a drop of quick-drying glue on the border between the *Cannula* and *Protrusion Adjustment Ring*.

## How to protect the *Cannula* when the *Input Protective Cap* does not stay inside?

1. Fill the interior of the *Cannula* with KWIK-CAST (WPI) to act as a cap. After removal of the dried sealant, clean the *Rod Lens* outer surface using a cotton swab lightly dipped in isopropyl alcohol.

## How do I stabilize the *Cannula-Microscope Body* connection?

1. Ensure that the Microscope Barrel is sufficiently loose (see *Microscope Surgical Procedure Application Note Section 2.1* or the *Instructional Video*).

## How to insert the *Microscope Body* into the *Cannula* with minimal force?

1. Ensure that the *Microscope Clamps* are sufficiently loose (screw in the barrel).
2. Ensure that the *Cannula* and *Microscope Body* are properly aligned.
3. Verify the *Cannula* installation instructions in *Section 2.1* of the *Microscope Surgical Procedure Application Note*.



### 7.0.3 Biology

Why can't I see any individual cells?

1. It generally takes 3 to 8 weeks before the area surrounding the tip of the *Rod Lens* has healed enough for sharp imaging of neurons.

## Support

### 8.1 Warranty

This product is under warranty for a period of 12 months. Contact Doric Lenses for return instructions. This warranty will not be applicable if the unit is damaged or needs to be repaired as a result of improper use or operation outside the conditions stated in this manual.

### 8.2 Contact us

For any questions or comments, do not hesitate to contact us by:

**Phone** 1-418-877-5600

**Email** [sales@doriclenses.com](mailto:sales@doriclenses.com)

The logo for Doric Lenses, featuring the word "doric" in a lowercase, sans-serif font. The letter 'o' is stylized with a white highlight on its upper-left curve, giving it a three-dimensional appearance.

© 2021 DORIC LENSES INC  
357 rue Franquet - Quebec, (Quebec)  
G1P 4N7, Canada  
Phone: 1-418-877-5600 - Fax: 1-418-877-1008  
[www.doriclenses.com](http://www.doriclenses.com)