



CAMERA DATA
(USB 3.0)

LED AND CAMERA CONTROL
(5X SMA)

EPHYS
(HDMI)

ROTATION CONTROL
(MINI USB2)

LASER INTERLOCK

POWER (12V)

ROTATION SENSOR

POWER INDICATOR

ON/OFF BUTTON

FOCUS
ADJUSTMENT

SAMPLE PORT (SMA)

FRONT VIEW

Ø 26.98

45.00

116.25

SIDE VIEW

1.25m

183.08

346.75

53.98

45.00

16 x Ø 3.45 ∇ 10.00
8-32 UNC ∇ 8.00

101.60

89.75

Ø 6.76

89.75

101.60

TOP VIEW

LIGHT SOURCE PORT
(M8 CONNECTOR)

FOR MORE DETAIL
SEE DRAWING:
#B340-2110-XXX

REAR VIEW

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF **Doric Lenses Inc.**
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING

doric

357 rue Franquet, Quebec, QC
G1P 4N7, Canada;
Tel: 418 877 5600;
doriclenses.com;

TITLE:

AHRJ-OE_24_RBFMC_G2

DWG. NO.

B330-1034-XXX-XX