



PATCHCORD
IMAGE FOR REPRESENTATION ONLY

CANNULA
SEE TABLE FOR PATTERN CONFIGURATION

CANNULA PATTERN CONFIGURATION TABLE							
LOCATION	TYPE	LENGTH	TIP	CANNULA	TYPE	LENGTH	TIP
C1-C10	Ø125 OR Ø250µm	0 - 10mm	FLT/Cxx/Axx/TPx.x	C11-C19	Ø125 OR Ø250µm	0 - 10mm	FLT/Cxx/Axx/TPx.x
C1				C11			
C2				C12			
C3				C13			
C4				C14			
C5				C15			
C6				C16			
C7				C17			
C8				C18			
C9				C19			
C10				-			

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF **Doric Lenses Inc.**..
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING

357 rue Franquet, Quebec, QC
G1P 4N7, Canada;
Tel: 418 877 5600;
doriclenses.com;

TITLE:

HD CANNULA PATTERN
CONFIGURATION DRAWING

DWG. NO.

N/A