


PATCHCORD CONFIGURATION TABLE			
FIBER TYPE AND LENGTH			
FIBER	PLASTIC 240/250-0.63 (ILLUMINATION ONLY) OR GLASS 200/220-0.37 (ILLUMINATION AND PHOTOMETRY RECORDING)		
LENGTH	STANDARD 2m		
LOCATION	PRESENCE	LOCATION	PRESENCE
P1-P10	YES OR NO	P11-P19	YES OR NO
P1		P11	
P2		P12	
P3		P13	
P4		P14	
P5		P15	
P6		P16	
P7		P17	
P8		P18	
P9		P19	
P10		-	

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF **Doric Lenses Inc.**..
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING



357 rue Franquet, Quebec, QC
G1P 4N7, Canada;
Tel: 418 877 5600;
doriclenses.com;

TITLE:

HD PATCHCORD PATTERN
CONFIGURATION DRAWING

DWG. NO.

N/A