



STRAIGHT PLOUT

MIRROR PINOUT

DATAMATE
12-WAY

120.0

22.6

10.0

18.3

4.4

11.6

7.5

2.0

doric

357 rue Franquet, Quebec, QC
G1P 4N7, Canada;
Tel: 418 877 5600;
doriclenses.com;

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF Doric Lenses Inc..
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING

TITLE:
**Harwin 12 to Omnetics
12 adapter**

DWG. NO.
B330-9006

8-1-2019