


8 7 6 5 4 3 2 1

D
C
B
A


D
C
B
A

USB TYPE B POWER SUPPLY JACK 2.5MM INTERLOCK BNC CONNECTORS (8X)






356

351



M8 CONNECTORS (4X)



8 7 6 5 4 3 2 1

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF **Doric Lenses Inc.**..
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING



357 rue Franquet, Quebec, QC
G1P 4N7, Canada;
Tel: 418 877 5600;
doriclenses.com;

TITLE:

LDMD_4

DWG. NO.

D480-2032

29/09/2017