




2	1	
FLT	DFL	C
Flat tip	Diffuser	Cone
		
Polished 0.3 μ m	Diffuser layer	

FIBER TERMINATION
SCALE 30 : 1

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF Doric Lenses Inc..
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING

	357 rue Franquet, Quebec, QC G1P 4N7, Canada; Tel: 418 877 5600; doriclenses.com;
TITLE:	OmFC_ZF1.25
DWG. NO.	B286-32xx
	06/06/2017