



ZIRCONIA

NICKEL PLATED BRASS
OR STAINLESS STEEL

HYTREL

$\phi 1.25$

5.0

11.9

$\phi 1.6$

<p>doric</p> <p>357 rue Franquet, Quebec, QC G1P 4N7, Canada; Tel: 418 877 5600; doriclenses.com;</p>	TITLE:
	<p>Generic MFP Termination for ZF1.25(RF) 1.6 mm-1100 and less</p>
	DWG. NO.
	<p>S212-0215</p>
<p>9-1-2025</p>	

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF **Doric Lenses Inc.**..
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING