


2	1	
FLT	DFL	C
Flat tip	Diffuser	Cone
Polished 0.3 µm	Diffuser layer	

FIBER TERMINATION
SCALE 50 : 1



DETAIL FIBER TERMINATION
SCALE 50 : 1

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF Doric Lenses Inc..
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING



357 rue Franquet, Quebec, QC
G1P 4N7, Canada;
Tel: 418 877 5600;
doriclenses.com;

TITLE:
OEC_ZF1.25

DWG. NO.
B290-xxxx

23-05-2017