

8 7 6 5 4 3 2 1

D

D

OPTICAL CONNECTOR

ELECTRICAL CONNECTOR

M3 CONNECTOR

142

Ø11


Ø6.35

C

C

B

B



A

A

8 7 6 5 4 3 2 1

Dimensions in mm  
For Information Purposes Only

**PROPRIETARY**  
THE INFORMATION CONTAINED IN THIS DRAWING  
IS THE SOLE PROPERTY OF **Doric Lenses Inc.**..  
ANY REPRODUCTION IN PART OR AS A WHOLE  
WITHOUT THE WRITTEN PERMISSION OF  
**Doric Lenses Inc.** IS PROHIBITED.

DO NOT SCALE DRAWING



357 rue Franquet, Quebec, QC  
G1P 4N7, Canada;  
Tel: 418 877 5600;  
doriclenses.com;

TITLE:  
**OEPH**

DWG. NO.  
B292-1xxx

10/11/2017